


Asepti Dry Wipes

A low lint, all-purpose towel as a companion item for Ecolab ready to use hard surface liquid cleaners which are suitable for all environmental cleaning and disinfectant applications



Features

- Asepti Dry Wipes are a perfect companion item for ECOLAB ready to use hard surface disinfectants.
- Asepti Dry Wipes are suitable for all environmental cleaning and disinfectant applications.
- Asepti Dry Wipes are non-sterile, single use, lint free dry wipes which are 35cm x 30cm in size and 55 grams per square.
- Asepti Dry Wipes are suitable for use in CSSD cleaning applications.
- The high absorbency of Asepti Dry Wipes are ideal for cleaning minor bodily fluids.
- Asepti Dry Wipes offer superior wet and dry strength and are latex free.

Product Code	Product Description	Pack Size	Quote
99001127, 30000372	Asepti Dry Wipes	8 x 100 wipes	

Safety Information

Do not mix with any other products. For further information and safety instructions see product label and Safety Data Sheet. ALWAYS READ THE LABEL, FOLLOW DIRECTIONS FOR USE.

Directions For Use

- Wear gloves and additional PPE where required.
- Apply to surface in a normal wiping motion.
- Dispose of wipe after use and always wash hands.
- Do not dispose of wipes in toilets.

Related Products

- [Asepti Active Liquid](#)
- [Asepti Active 24/7](#)
- [Actichlor Plus](#)

Applications

- Infection Control
- Residential Aged Care
- Healthcare Facilities
- CSD
- Veterinary
- Wards
- Theatre
- Dental
- ICU
- Day Surgeries
- Patient Rooms
- Bathrooms

Resources

- Video
- Brochure