

Belimed SmartHub

Belimed SmartHub is your gateway to a smart, innovative CSD.



Features

- Belimed SmartHub is an on-premise software that ensures all relevant data from your Belimed equipment is captured in real time to be further processed by either an Instrument Tracking System or via an easy, readable PDF.
- Belimed reprocessing equipment produce a huge amount of data points per cycle, including items like program selection, chemical dosing quantities, temperature, water volume, operating times and water conductivity. This information is captured and processed through Belimed Smarthub, pushing through to a dashboard system for you to analyse and make informed decisions.
- Belimed Smarthub is your gateway to the cloud and digitalising your CSD. Digitalising leads to a transparent CSD, enabling a smart CSD, allowing the Belimed machines to directly communicate with the end user anytime and anywhere throughout the facility.
- Belimed Smarthub operates as a data turntable that is used for; Receiving machine and cycle data from the connected devices, Temporarily storing the data within a database, Transforming the data into a required format, Providing paperless storage and retrieval of reprocessing documentation, Exporting the data to a file system on a local server and As an enabler for the digital CSD.
- Being part of a critical infrastructure, cybersecurity is of Belimed’s utmost priority. Belimed developed the connect software products inspired by the standard IEC 82304-1 (health software) and in development underwent in-depth extensive testing according to the agile testing quadrants of Crispin/Marick. While Belimed SmartHub does not process any patient or operator data, Belimed provide the option of encrypting communication within the hospital as well as the outside world.

Product Code	Product Description	Pack Size	Quote
BMA-1053485, SER99985	Belimed SmartHub Site Licence	ea	
BMA-1053465, SER90496	Belimed SmartHub Connect	ea	
BMA-1053480, SER90497	Belimed SmartHub Dashboard	ea	
BMA-1053481, SER90535	Belimed SmartHub Machine Connections	ea	
BMA-1053487, SER90490	Belimed SmartHub Connect, Upgrade ICS	ea	

Product Code	Product Description	Pack Size	Quote
BMA-1053489, SER90491	Belimed SmartHub Connect, Upgrade SH1	ea	
BMA-1056617, SER90492	Belimed SmartHub Multisite, Upgrade ICS	ea	
BMA-1056619, SER90493	Belimed SmartHub Dashboard, Upgrade ICS	ea	
BMA-1056620, SER90494	Belimed SmartHub Multisite, Upgrade SH1	ea	
BMA-1056621, SER90495	Belimed SmartHub Dashboard, Upgrade SH1	ea	
BMA-1064346, SER90498	Belimed SmartHub Orbit Annual Licence	ea	
BMA-1064505, SER90534	Belimed SmartHub Orbit - Advanced Subscription	ea	
DIG-INSTALL- GENMED, EXP01569	Deployment & Testing of SmartHub Software (remote)	ea	

Relevant Product

- [WD290IQ Washer Disinfector](#)
- [WD290 Washer Disinfector](#)
- [WD200 Washer Disinfector](#)
- [WD390 Index Washer](#)
- [WD750 Cart Washer](#)
- [CS750 Clean Station](#)
- [MST-V Steam Steriliser](#)
- [MST-H Steam Steriliser](#)
- [DC2200 iX Drying Cabinet](#)
- [UW565 Ultrasonic](#)

Applications

- CSD
- Theatre
- Endoscopy
- Healthcare Facilities
- Dental
- Life Science

Resources

- video
- video
- video
- video
- video
- Brochure