


## Belimed SmartHub

Belimed SmartHub is your gateway to a smart, innovative CSD.



## Features

- Belimed SmartHub is an on-premise software that ensures all relevant data from your Belimed equipment is captured in real time to be further processed by either an Instrument Tracking System or via an easy, readable PDF.
- Belimed reprocessing equipment produce a huge amount of data points per cycle, including items like program selection, chemical dosing quantities, temperature, water volume, operating times and water conductivity. This information is captured and processed through Belimed Smarthub, pushing through to a dashboard system for you to analyse and make informed decisions.
- Belimed Smarthub is your gateway to the cloud and digitalising your CSD. Digitalising leads to a transparent CSD, enabling a smart CSD, allowing the Belimed machines to directly communicate with the end user anytime and anywhere throughout the facility.
- Belimed Smarthub operates as a data turntable that is used for; Receiving machine and cycle data from the connected devices, Temporarily storing the data within a database, Transforming the data into a required format, Providing paperless storage and retrieval of reprocessing documentation, Exporting the data to a file system on a local server and As an enabler for the digital CSD.
- Being part of a critical infrastructure, cybersecurity is of Belimed's utmost priority. Belimed developed the connect software products inspired by the standard IEC 82304-1 (health software) and in development underwent in-depth extensive testing according to the agile testing quadrants of Crispin/Marick. While Belimed SmartHub does not process any patient or operator data, Belimed provide the option of encrypting communication within the hospital as well as the outside world.

Product Code	Product Description	Pack Size	Quote
1040274	Belimed SmartHub	ea	

## Relevant Product

- [WD290IQ Washer Disinfector](#)
- [WD290 Washer Disinfector](#)
- [WD200 Washer Disinfector](#)
- [WD390 Index Washer](#)
- [WD750 Cart Washer](#)
- [CS750 Clean Station](#)

- [MST-V Steam Steriliser](#)
- [MST-H Steam Steriliser](#)
- [DC2200 iX Drying Cabinet](#)
- [UW565 Ultrasonic](#)

## Applications

---

- CSD
- Theatre
- Endoscopy
- Healthcare Facilities
- Dental
- Life Science

## Resources

---

- video
- video
- video
- video
- video
- Brochure