

Belimed WD750

Large-scale cleaning, disinfection and drying system for safe and efficient cleaning of large laboratory and pharmaceutical equipment

ARTG R 228756



Features

- Belimed WD750 large scale washer disinfectors sets new standards in safety, quality and economic efficiency.
- Belimed WD750 is available in 3 sizes depending on your CSD requirements; small, medium or large chamber.
- Belimed WD750 large scale washer disinfectors support your workflow through high quality components and processing, integrated, innovative and patented features along with the machines economic design and performance.
- Belimed WD750 is the most compact large-scale washer in its class, incorporating easy service access and integrated equipment/media connections accommodating media connections and detergents all in one area.
- Belimed WD750's glass, horizontal sliding door, provides clear visibility into wash chamber; 1 or 2 door versions available.
- The patented technology, LED status display, provides easy use and progress information on cycle progress.
- Belimed WD750 large scale washer disinfectors access a versatile range of racks, specifically suited to the items to be washed.

Product Code	Product Description	Pack Size	Quote
BMC-100-2926, 10119373	Belimed WD750, Small, 2 Door	ea	
BMC-100-2928, 10114420	Belimed WD750, Medium, 2 Door	ea	
BMC-100-2930, 10114421	Belimed WD750, Large, 2 Door	ea	

Technical Data

Dimensions (mm)	WD750S	WD750M	WD750L
Product Outer (H x W x D)	2450 x 2300 x 2000	2450 x 2300 x 2860	2450 x 2300 x 3600
Product Chamber (H x W x D)	2000 x 820 x 1350	2000 x 820 x 2300	2000 x 820 x 2950
Chamber Volume (L)	2214	3772	4838

Load Capacities			
Sterile Containers with Lid (300 x 600 x 300mm)	20	27	40
Sterile Cages for Utensils (340 x 570 x 185mm)	16	16	32
Surgical Instruments DIN Sieve Trays	28	44	56
Threatre Dishes Ø 320mm	28	44	56
Threatre Shoes up to size 42	180	270	360
Theatre Shoes size 43-46	108	162	216
Carts	1	2	2
Threatre Tables	0	1	1

Relevant Product

- [Asepti RHW](#)
- [Asepti NewGen](#)
- [Asepti Alcohol Rinse](#)
- [Asepti Pro-Dry](#)
- [Belimed Washer Disinfector](#)
- [Belimed Large Scale Washer](#)
- [Terragene Thermal Disinfection Control](#)

Applications

- Life Science

Resources

- Video
- Video
- Video
- Brochure