

## Belimed WD750

Large-scale cleaning, disinfection and drying system for safe and efficient cleaning of medical and non-medical devices

ARTG 228756



## Features

- Belimed WD750 large scale washer disinfectors sets new standards in safety, quality and economic efficiency.
- Belimed WD750 is available in 3 sizes depending on your CSD requirements; small, medium or large chamber.
- Belimed WD750 large scale washer disinfectors support your workflow through high quality components and processing, integrated, innovative and patented features along with the machines economic design and performance.
- Belimed WD750 is the most compact large-scale washer in its class, incorporating easy service access and integrated equipment/media connections accommodating media connections and detergents all in one area.
- Belimed WD750's glass, horizontal sliding door, provides clear visibility into wash chamber; 1 or 2 door versions available.
- The patented technology, LED status display, provides easy use and progress information on cycle progress.
- Belimed WD750 large scale washer disinfectors access a versatile range of racks, specifically suited to the items to be washed.

| Product Code              | Product Description           | Pack Size | Quote                                                                                 |
|---------------------------|-------------------------------|-----------|---------------------------------------------------------------------------------------|
| BMC-100-2926,<br>10119373 | Belimed WD750, Small, 2 Door  | ea        |  |
| BMC-100-2928,<br>10114420 | Belimed WD750, Medium, 2 Door | ea        |  |
| BMC-100-2930,<br>10114421 | Belimed WD750, Large, 2 Door  | ea        |  |

## Technical Data

| Dimensions (mm)             | WD750S             | WD750M             | WD750L             |
|-----------------------------|--------------------|--------------------|--------------------|
| Product Outer (H x W x D)   | 2450 x 2300 x 2000 | 2450 x 2300 x 2860 | 2450 x 2300 x 3600 |
| Product Chamber (H x W x D) | 2000 x 820 x 1350  | 2000 x 820 x 2300  | 2000 x 820 x 2950  |
| Chamber Volume (L)          | 2214               | 3772               | 4838               |

| Load Capacities                                 |     |     |     |
|-------------------------------------------------|-----|-----|-----|
| Sterile Containers with Lid (300 x 600 x 300mm) | 20  | 27  | 40  |
| Sterile Cages for Utensils (340 x 570 x 185mm)  | 16  | 16  | 32  |
| Surgical Instruments DIN Sieve Trays            | 28  | 44  | 56  |
| Theatre Dishes Ø 320mm                          | 28  | 44  | 56  |
| Theatre Shoes up to size 42                     | 180 | 270 | 360 |
| Theatre Shoes size 43-46                        | 108 | 162 | 216 |
| Carts                                           | 1   | 2   | 2   |
| Theatre Tables                                  | 0   | 1   | 1   |

## Relevant Product

- [WD290IQ Washer Disinfector](#)
- [WD290 Washer Disinfector](#)
- [CS750 Washer Disinfector](#)
- [WD7015 Washer Disinfector](#)
- [MST-V Steam Steriliser](#)
- [MST-H Steam Steriliser](#)
- [Belimed Smarthub](#)
- [CSD Accessories](#)
- [Asepti Multizyme](#)
- [Asepti Neutrazyme](#)
- [Asepti Synergy WD](#)
- [Central Service Program](#)

## Applications

- CSD
- Healthcare Facilities
- Large Dental Practices
- Life Science

## Resources

- Video
- Video
- Video
- Video
- Brochure