


Dermados Dispenser 1L

Dispenser to suit 1 L pump



Features

- Dermados dispensers offer a reliable and economical dispensing solution.
- Touch free elbow operation.
- Dermados dispensers provide a secure application of alcohol-based hygienic and surgical rubs in all areas.
- Easy installation with and without drilling option.
- Volume adjustable from 0.5 mL to 1.5 mL.
- Dispenser is easy to clean and to disinfect using commercial cleaning and disinfecting agents and is dishwasher safe up to 70°C.
- Manufactured using impact resistant ABS plastic.
- Corrosion resistant against commercially available cleaning and disinfecting products.
- Low maintenance solution.
- Lockable for peace of mind, especially for dispensers located in public areas.

Product Code	Product Description	Pack Size	Quote
9900-1002, 10056798	Dermados Dispenser to suit 1 L Pump	ea	

Relevant Product

- [Skinman 90, 1 L](#)

Resources

- Video
- Brochure