


## Distal Tip Protectors

Flexible Endoscope Tip Protectors

ARTG 227568



### Features

- EndoSelect flexible endoscope distal tip protectors are the ideal solution for protecting the scopes tip during transport.
- EndoSelect Distal Tip Protectors are made from semi-rigid expandable mesh to protect the tip/lense of the instrument from front and side impact.
- They are single use and fit the common Gastrosopes, Colonoscopes, Duodenoscopes brands on the market.
- EndoSelect Distal Tip Protectors allow moisture to escape, eliminating potential bacterial growth.
- EndoSelect Distal Tip Protectors are designed not to fall off and are non-abrasive so as not to scratch any part of the instrument or lens.
- EndoSelect Distal Tip Red: Small for instruments approximately 2-10mm diameter
- EndoSelect Distal Tip Blue: Large for instruments approximately 7-15mm diameter

Product Code	Product Description	Pack Size	Quote
GM-TIP550100, 10116924	EndoSelect Small Red Tip Protectors	100 Pk	
GM-TIP650100, 10116925	EndoSelect Large Blue Tip Protectors	100 Pk	

### Safety Information

Do not mix with any other products. Further information, safety instructions and health warnings, see product label and Safety Data Sheet.  
 ALWAYS READ THE LABEL. FOLLOW DIRECTIONS OF USE.

### Relevant Product

- [Soluscope Sprint](#)
- [Soluscope Sprint PT](#)
- [Soluscope DSC8000](#)
- [EndoSelect Cleaning Sponges](#)
- [EndoSelect Buttons & Valves](#)
- [EndoSelect Patient Pad](#)
- [EndoSelect Quad Trap](#)
- [Soluscope Pro Specimen Trap](#)

- [Asepti Epizyme Ultra Wipes](#)
- [Cleaning Bowls](#)
- [EndoSelect Cleaning Brushes](#)
- [EndoSelect Transport Pad](#)
- [EndoSelect Bite Blocks](#)
- [EndoSelect Distal Tip Protectors](#)
- [EndoSelect Gauze](#)

## Applications

---

- Endoscopy
- Healthcare Facilities

## Resources

---

- Video
- Brochure