






## Smartline SmartSink





















Electric Height Adjustable Washing Stations improve productivity












## Features

- Smartline SmartSink offers your reprocessing department an electric height adjustable washing station which improves productivity and reduces the risk of staff injury.
- Smartline SmartSinks are manufactured in Australia and constructed out of 304 grade stainless steel. The Smartline range of SmartSinks offer a complete, customisable washing station solution for your departments requirements.
- Smartline SmartSinks feature a 300mm high tapered splashback and tapered benchtop with wetdeck perimeter, offering smooth surfaces for easy cleaning for infection prevention.
- SmartSinks include easy-to-use electric height adjustment with anti-collision technology from 850mm to 1080mm off the ground. Controls are positioned under the sink top offering the user full control and easy access.
- Smartline SmartSinks offer so many customisable options; multiple bench surface, two different or customisable bowl depth, rectangular, circular, single and/or double sink options, sink inserts, shelves, trays and racks and tapware to name a few.
- All Smartline SmartSinks are manufactured with adjustable bolt down feet.

Product Code	Product Description	Pack Size	Quote
11.SS-135-SC	Smartline 1350mm SmartSink Single Bowl - Circular	ea	
11.SS-175-SC	Smartline 1750mm SmartSink Single Bowl - Circular	ea	
11.SS-195-SC	Smartline 1950mm SmartSink Single Bowl - Circular	ea	
11.SS-220-SC	Smartline 2200mm SmartSink Single Bowl - Circular	ea	

Product Code	Product Description	Pack Size	Quote
11.SS-175-DC	Smartline 1750mm SmartSink Double Bowl - Circular	ea	
11.SS-195-DC	Smartline 1950mm SmartSink Double Bowl - Circular	ea	
11.SS-220-DC	Smartline 2200mm SmartSink Double Bowl - Circular	ea	
11.SS-260-DC	Smartline 2600mm SmartSink Double Bowl - Circular	ea	
11.SS-320-DC	Smartline 3200mm SmartSink Double Bowl - Circular	ea	
11.SS-360-DC	Smartline 3600mm SmartSink Double Bowl - Circular	ea	
11.SS-UT13560	Smartline Under Sink Tray - suit 1350x600 Sink Units	ea	
11.SS-UT13570	Smartline Under Sink Tray - suit 1350x700 Sink Units	ea	
11.SS-UT17560	Smartline Under Sink Tray - suit 1750x600 Sink Units	ea	
11.SS-UT17570	Smartline Under Sink Tray - suit 1750x700 Sink Units	ea	
11.SS-UT19570	Smartline Under Sink Tray - suit 1950x700 Sink Units	ea	
11.SS-UT22060	Smartline Under Sink Tray - suit 2200x600 Sink Units	ea	
11.SS-UT22070	Smartline Under Sink Tray - suit 2200x700 Sink Units	ea	
11.SS-UT26060	Smartline Under Sink Tray - suit 2600x600 Sink Units	ea	
11.SS-UT26070	Smartline Under Sink Tray - suit 2600x700 Sink Units	ea	
11.SBS2050	Smartline Smartsink Splash Back Shelf - 20x50cm	ea	
11.SBS2080	Smartline Smartsink Splash Back Shelf - 20x80cm	ea	
11.SBT2050	Smartline SmartSink Splashback Tray 20x50cm	ea	
11.SBT2080	Smartline SmartSink Splashback Tray 20x80cm	ea	
11.SS-CHEMSHLF	Smartline SmartSink Chemical Bottle Shelf	ea	

Product Code	Product Description	Pack Size	Quote
11.SS-SBCM	Smartline SmartSink Clamp Mount for Accessories	ea	
51.TRYRCK-1X5	Smartline SmartSink Stainless Steel Tray Rack for Autoplas Trays	ea	
51.TRDIN-1X5	Smartline SmartSink Stainless Steel Tray Rack for DIN Trays	ea	
11.SS-GF152510-SS-SH	Smartline SmartSink Pre-Rinse Unit, 10 inch spout	ea	
11.TAP-CL10010	Smartline SmartSink Lever Mixing Set	ea	
11.TAP-CL15035	Smartline SmartSink Lever Sink Set	ea	
11.TAP-JETF3000W	Smartline SmartSink Pre-Rinse Unit	ea	
47.DINA-JET	Smartline SmartSink Cleaning Jet Pistol with Standard Set	ea	
47.HELTA-JET	Smartline SmartSink Cleaning Jet Pistol with Standard Set, Stainless Steel	ea	

## Relevant Product

- [Asepti Epizyme Ultra](#)
- [Asepti Multizyme](#)
- [Asepti Neutrazyme](#)
- [Asepti Rapidzyme PR](#)
- [Asepti Release Plus](#)
- [Asepti Synergy 5](#)
- [Asepti Synergy WD](#)
- [Asepti Wash](#)
- [Optipro Multi Enzyme](#)
- [Water Testing, Table 7.3](#)
- [EndoSelect Cleaning Brushes](#)
- [EndoSelect Cleaning Sponges](#)
- [Transportation Trolleys](#)
- [Central Service Program](#)

## Applications

- Endoscopy
- Central Sterile Department
- Healthcare Facilities

## Resources

- video
- Brochure